

M30 CytoDeath™ ELISA

Apoptosis detection in 2D and 3D cell cultures

Catalog Prod. No. 10900

In USA, Canada and Japan: For research and laboratory use only. Not for human or diagnostic use.

General Information

Analyte:	Soluble caspase-cleaved fragments of the intermediate filament protein keratin 18 (K18) containing the M30 neo-epitope (K18Asp396-NE). The assay detects human, monkey and bovine K18 fragments.
Intended Use:	<p>Quantitative measurement of the apoptotic cell death biomarker K18Asp396-NE in cell culture experiments. Can be used for cell lysates and/or culture supernatants. The assay only detects apoptosis in cells of epithelial origin that express K18. Cells should be of human, monkey or bovine origin.</p> <p>To be used to determine accumulation of caspase-cleaved K18 in cell cultures, providing an integrative measure of apoptosis. The K18Asp396 neo-epitope is formed by caspase-3, -7 or -9 activation.</p>
Samples:	Cell lysates or culture supernatants. <i>Not suitable</i> for serum or plasma samples.
Sample Volume:	2 × 25 µl (duplicate samples).
Sample Stability:	Fresh samples are stable for up to two days at 2–8 °C, for at least 9 months at -20 °C; and for at least two years when stored at -80 °C.
Number of Tests:	96 determinations: 4 Standards and 44 samples in duplicates.
Reagent Storage:	2 – 8 °C. Do not freeze!
Assay Time:	260 min (approx.).
References:	<ul style="list-style-type: none">■ Fayad W, et al. (2009). <i>Identification of a novel topoisomerase inhibitor effective in cells overexpressing drug efflux transporters</i>. PLoS One 4(10):e7238.■ Hernlund E, et al. (2009). <i>Ovarian carcinoma cells with low levels of b-F1-ATPase are sensitive to combined platinum and 2-deoxy-D-glucose treatment</i>. Mol Cancer Ther 8(7).■ Herrmann R, et al. (2008). <i>Screening for Compounds that Induce Apoptosis of Cancer Cells Grown as Multicellular Spheroids</i>. J Biomol Screen. 13(1):1 – 8.■ Lakshmikanthan V, et al. (2006). <i>SAHA-sensitized prostate cancer cells to TNFalpha-related apoptosis-inducing ligand (TRAIL): mechanisms leading to synergistic apoptosis</i>. Int J Cancer 119:221 – 8.■ Erdal H, et al. (2005). <i>Induction of lysosomal membrane permeabilization by compounds that activate p53-independent apoptosis</i>. Proc. Natl. Acad. Sci. USA 102, 192 – 197.■ Schutte B, et al. (2004). <i>Keratin 8/18 breakdown and reorganization during apoptosis</i>. Exp Cell Res. 297, 11 – 26.■ Kramer G, et al. (2004). <i>Differentiation between Cell Death Modes using Measurements of Different Soluble Forms of Extracellular Cytokeratin 18</i>. Cancer Research 64, 1751 – 1756.■ Hägg, M. et al. (2002). <i>A novel high-through-put assay for screening of pro-apoptotic drugs</i>. Invest. New Drugs, 20: 253-259.■ Leers MP, et al. (1999). <i>Immunocytochemical detection and mapping of a Cytokeratin 18 neo-epitope exposed during early apoptosis</i>. J Pathol 187, 567-572.

Performance Characteristics

Calibration:	The Units measured by the M30 CytoDeath™ ELISA are defined against a synthetic peptide containing the M30 and M6 epitopes. 1 U/L = 1.24 pM.
Working Range:	250 – 3 000 U/L.
Detection Limit:	60 U/L, Standard Z (0 U/L) + 3 S.D.
Reproducibility:	Intra-Assay (WA) Precision: CV < 7 % for values > 250 U/L. Inter-Assay (BA) Precision: CV < 10 % for values > 250 U/L.
Spike Recovery:	80–120 %
Linearity/Dilution:	80–120 %
Hook Effect:	No high dose “hook effect” occurs before 26 000 U/L which is well above concentrations of K18Asp396-NE (M30)-reactive material observed in cell culture samples.

Reagents

Coated Microstrips:	One Microplate, 96 dry wells (12 strips × 8 wells). The wells are coated with mouse monoclonal K18 anti-body M6.
HRP Conjugate:	Concentrate. One vial containing 0.4 mL mouse monoclonal M30 antibody (anti-K18Asp396-NE) conjugated to horseradish peroxidase (HRP).
Conjugate Dilution Buffer:	One vial containing 11 mL of phosphate buffer with protein stabilizers.
Standards:	The Standards are Standard Zero (0 U/L), Low (250 U/L), Medium (1 000 U/L) and High (3 000 U/L).
TMB Substrate:	One vial containing 22 mL of TMB (3,3',5,5'-Tetramethylbenzidine) Solution.
Stop Solution:	One vial containing 6 mL of 1.0 M sulfuric acid.
Wash Tablet:	One tablet for preparation of 500 mL Wash Solution.

Assay Advantages

Easy to use:	All reagents are ready-to-use except for the M30-HRP Conjugate concentrate and the Wash Tablet. Only few steps and a limited hands-on time are required. Results are recorded as absorbance at 450nm using standard 96 well plate readers and are quantified against a standard.
Unique application range:	Coated Microstrips: One Microplate, 96 dry wells (12 strips × 8 wells). The wells are coated with mouse monoclonal K18 anti-body M6. antibody Detecting specifically apoptosis, M30 CytoDeath™ ELISA does not reflect cellular metabolic activity often measured by viability tests as this may also be compromised during drug-induced cell cycle arrest, senescence or differentiation. The assay shows a unique application to selectively detect apoptotic K18-positive target cells (e.g. carcinoma) in co-cultures with cells that are negative for K18 such as fibroblasts, endothelial cells and macrophages (e.g. stroma) or 3D tissue/tumor cultures and multi-cellular tumor spheres.

PEVIVA Products from VLVbio

M30 Apoptosense® ELISA Prod. no. 10011	M65® ELISA Prod. no. 10020	M5 Keratin 18 Prod. no. 10600	M30 CytoDEATH™ Unconjugated Prod. No. 10700 Biotin Prod. No. 10750 Fluorescein Prod. No. 10800 Orange Prod. No. 10850
M30 CytoDeath™ ELISA Prod. no. 10900	M65 EpiDeath® ELISA Prod. no. 10040	M6 Keratin 18 Prod. no. 10650	

For further information, please visit www.peviva.com - order online at www.shop.peviva.co.uk

VLVbio 

Vivalavida AB
172 66 Sundbyberg
Sweden
info@vlvbio.com
www.vlvbio.com



Intertek

bioaxxess


BIOAXXESS
Basepoint Business Centre
GL20 8SD Tewkesbury
United Kingdom

Your Partner in UK

Phone +44 (0) 1684 851 270
Fax +44 (0) 1684 851 272
Email info@bioaxxess.com
Web www.bioaxxess.co.uk